
## Revisions

SO 0191-5

Rev from	Rev to	ECO #
0603	0105	3277-05

**Notes:**

1. BD Cat. Number 223001, 223011, 223021, 223031, 223041, 223051, 223061
2. Blank (Sheet) Size : Length: 9"      Width: 36"  
 Number of Pages: 18      Number of Sheets: 1  
 Page Size: Length 9"      Width 4"      Final Folded Size: 2.25" x 4"
3. Style (see illustrations below): #4

4. See Specification Control Number S1471 for Material Information
5. Ink Colors: Printed two sides  Yes     No  
 No. of Colors: 1      PMS #2755 Blue
6. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	<b>Becton, Dickinson and Company</b> 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number:		Category and Description Package Insert Listeria Antigens and Antisera	Sheet: 1 of 19	A
S1471			Scale: 1:1	

# BD Difco™ Listeria Antisera and Antigens

 S1471  
2005/01

**English:** pages 1 - 4      **Italiano:** pagine 11 - 13  
**Français :** pages 4 - 7      **Español:** páginas 14 - 17  
**Deutsch:** Seiten 7 - 10

See symbol glossary at end of insert. / Viz popis symbolů na konci příbalového letáku. / Se symbolglossaret i slutningen af indlægssedlen. / Zie lijst met symbolen aan het einde van de bijsluiter. / Vaadake sümbolite seletust infolehe lõpus. / Katso pakkauselosteen lopussa olevaa kuvamerkkien sanastoa. / Voir le glossaire des symboles à la fin de la notice. / Siehe Symbol-Erklärungen am Ende der Packungsbeilage. / Δείτε το γλωσσάριο των συμβόλων στο τέλος του ένθετου. / A jelmagyarázat a használati utasítás végén található. / Vedere il glossario dei simboli alla fine del foglio illustrativo. / Zr. informacino lapelio pabaigoje pateikiama simbolių glosarijų. / Se i symbolforklaringen på slutten av produktvedlegget. / Zobacz objaśnienie symboli na końcu ulotki. / Consulte o glossário de símbolos no fim do folheto informativo. / Pozri slovník symbolov na konci letáka. / Consulte el glosario de símbolos al final del prospecto. / Se symbolförteckningen vid slutet av bipacksedeln.

Pokyny vám poskytnú miestni zástupce spoločnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselőjétől. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o seu representante local da BD para obter instruções. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar.

## INTENDED USE

**Difco™** Listeria O Antisera Types 1, 4 and Poly are used for identifying *Listeria monocytogenes* in the slide and tube agglutination tests.

**Difco** Listeria O Antigens Types 1 and 4 (Slide and Tube) are used as positive controls in the slide and tube agglutination tests, respectively.

## SUMMARY AND EXPLANATION

First described in 1926 by Murray, Webb and Swann,<sup>1</sup> *Listeria monocytogenes* is a widespread problem in public health and the food industries. This organism can cause human illness and death, particularly in immunocompromised individuals and pregnant women.<sup>2</sup> The first reported foodborne outbreak of listeriosis was in 1985<sup>3</sup> and, since then, microbiological and epidemiological evidence from both sporadic and epidemic cases of listeriosis has shown that the principal route of transmission is via the consumption of foodstuffs contaminated with *L. monocytogenes*.<sup>4</sup>

The organism has been isolated from commercial dairy and other food processing plants and is ubiquitous in nature, being present in a wide-range of unprocessed foods and in soil, sewage, silage and river water.<sup>5</sup>

*Listeria* spp. grow over a pH range of 4.4-9.6 and survive in food products with pH levels outside these parameters.<sup>6</sup> *Listeria* spp. are microaerophilic, gram-positive, asporogenous, non-encapsulated, non-branching, regular, short, motile rods. Motility is most pronounced at 20°C.

Strains of *Listeria* spp. are divided into serotypes based on cellular (O) and flagellar (H) antigens.<sup>7</sup> Thirteen serotypes of *L. monocytogenes* are known. Most human disease is caused by serotypes 1/2a, 1/2b and 4b.<sup>8</sup>

## PRINCIPLES OF THE PROCEDURE

Identification of *L. monocytogenes* includes both biochemical and serological confirmation.

Serological confirmation requires that the microorganism (antigen) react with its corresponding antibody. This *in vitro* reaction produces macroscopic clumping called agglutination. The desired homologous reaction is rapid, does not dissociate (high avidity) and binds strongly (high affinity).

Because a microorganism (antigen) may agglutinate with an antibody produced in response to another species, heterologous reactions are possible. These are characterized as weak in strength or slow in formation. Such unexpected and perhaps unpredictable reactions may lead to some confusion in serological identification. A positive homologous agglutination reaction should support the morphological and biochemical identification of the microorganism.

Agglutination of the somatic antigen in the slide test appears as a firm granular clumping. Homologous reactions occur rapidly and are strong (3+). Heterologous reactions form slowly and are weak.

The agglutination of the somatic antigen in the tube test appears as a loose flocculation that can easily be resuspended. Homologous reactions using Listeria O Antisera should exceed a titer of 2+ at 1:320.

## REAGENTS

**Difco** Listeria O Antisera Types 1, 4 and Poly are lyophilized, polyclonal rabbit antisera containing approximately 0.04% Thimerosal as a preservative. Listeria O Antisera Types 1 and 4 are specific for the respective serotypes of *L. monocytogenes* while Listeria O Antiserum Poly contains agglutinins for *L. monocytogenes* serotypes 1 and 4.

**Difco** Listeria O Antigens Types 1 and 4 (Tube) and (Slide) are suspensions of appropriate *L. monocytogenes* serotypes containing 0.3% formaldehyde as a preservative. When used according to the suggested procedure, the reagents will yield the following:

Reagent	Vial	Number of tests
Listeria O Antiserum	1 mL	10 tube tests or 400 slide tests
Listeria O Antigen (Slide)	5 mL	100 slide tests
Listeria O Antigen (Tube)	25 mL	5 tube tests

#### Warnings and Precautions

For *in vitro* Diagnostic Use.

The Packaging of This Product Contains Dry Natural Rubber.

Observe aseptic techniques and established precautions against microbiological hazards throughout all procedures. After use, specimens, containers, slides, tubes and other contaminated material must be sterilized by autoclaving. Directions for use should be followed carefully.

**Storage:** Store lyophilized and rehydrated **Difco** Listeria O Antisera at 2-8°C.

Store **Difco** Listeria O Antigen (Slide) and (Tube) at 2-8°C.

Prolonged exposure of reagents to temperatures other than those specified is detrimental to the products.

**Product Deterioration:** The expiration date applies to the product in its intact container when stored as directed. Do not use if the product is caked, discolored or shows other signs of deterioration. Examine antigen vials for agglutination before use. Suspensions with agglutination are not usable and should be discarded. Discard any antiserum that is cloudy or has a precipitate after rehydration or storage.

#### SPECIMEN COLLECTION AND PREPARATION

**Clinical specimens:** *Listeria* can be recovered on selective differential media such as McBride Listeria Agar, Oxford Agar, Modified Oxford Agar, LPM Agar or PALCAM Medium. For specific recommendations on isolation of *Listeria* from clinical specimens, consult appropriate references.<sup>8-10</sup> Determine that a pure culture of the microorganism has been obtained and that biochemical test reactions are consistent with the identification of the organism as *L. monocytogenes*. After these criteria are met, serological identification can be performed.

**Food or dairy samples:** *Listeria* can be recovered when samples are processed to recover injured microorganisms and prevent overgrowth of competing microorganisms. Consult appropriate references for recommended procedures for the isolation of *Listeria* from foods.<sup>6,11,12</sup> Having followed an established protocol, isolate a pure culture of the microorganism and confirm that biochemical test reactions are consistent with the identification of the organism as *L. monocytogenes*. After these criteria are met, serological identification can be performed.

#### PROCEDURE

**Materials Provided:** **Difco** Listeria O Antiserum Type 1, **Difco** Listeria O Antiserum Type 4, **Difco** Listeria O Antiserum Poly, **Difco** Listeria O Antigen Type 1 (Slide), **Difco** Listeria O Antigen Type 1 (Tube), **Difco** Listeria O Antigen Type 4 (Slide), **Difco** Listeria O Antigen Type 4 (Tube).

#### Materials Required But Not Provided:

**Slide Test:** **Difco** FA Buffer (Dried), agglutination slides, applicator sticks, water bath (80-100°C), droppers.

**Tube Test:** **Difco** FA Buffer (Dried), McFarland turbidity standard no. 3, culture tubes 12 x 75 mm and rack, serological pipettes, 1 mL, water bath 50°C, refrigerator (2-8°C), formaldehyde.

**Reagent Preparation:** Equilibrate all materials to room temperature before performing the tests. Ensure that all glassware and pipettes are clean and free of residues such as detergents.

**Listeria O Antisera:** To rehydrate, add 1 mL sterile purified water to each vial. Rotate gently to dissolve contents completely.

**Listeria O Antigens** (Slide) and (Tube) are ready to use.

#### Test Procedure

##### Slide Test

1. Rehydrate **Difco** FA Buffer as per label directions.
2. Suspend a loopful of the test isolate growth from a solid noninhibitory agar medium in approximately 5 mL of **Difco** FA Buffer.
3. Heat the organism suspension at 80-100°C (in a water bath) for 1 h.
4. Centrifuge the suspension at 2,000-5,000 RCF for 15-20 min and remove the bulk of the supernatant fluid.
5. Resuspend the organism in the remaining portion of supernatant fluid.
6. Prepare a 1:20 dilution of the desired **Difco** Listeria Antiserum in 0.85% NaCl solution.
7. On an agglutination slide, dispense 1 drop of the desired diluted **Difco** Listeria Antiserum onto each of two separate areas. Dispense 1 drop of **Difco** FA Buffer onto a third area of the same slide. The first drop of antiserum will be used for the test isolate and the second for the positive control. The drop of **Difco** FA Buffer will be used for the negative control.
8. Dispense 1 drop of organism suspension from step 5 to the first drop of antiserum.
9. **Positive control:** Dispense 1 drop of homologous **Difco** Listeria O Antigen (Slide) to the second drop of antiserum.

10. **Negative control:** Dispense 1 drop of organism suspension from step 5 to the drop of **Difco** FA Buffer on the agglutination slide.
11. Rotate the slide for 1-2 min and read for agglutination.

#### Tube Test

1. Suspend a loopful of the test organism from a solid agar medium in formalinized **Difco** FA Buffer (0.3 mL formaldehyde per 300 mL **Difco** FA Buffer). Adjust to a density approximating that of a McFarland turbidity standard no. 3.
2. Prepare a row of 9 culture tubes (12 x 75 mm) for each organism suspension to be tested, including one row for the positive control.
3. Dispense 0.9 mL formalinized **Difco** FA Buffer to the first tube in each row and 0.5 mL to the remaining tubes.
4. Using a 1 mL serological pipette, add 0.1 mL of the desired **Difco** Listeria Antiserum to tube 1 in each row and mix thoroughly. Transfer 0.5 mL from tube 1 to tube 2 and mix thoroughly. In a like manner, continue transferring 0.5 mL through tube 8, discarding 0.5 mL from tube 8 after mixing. Tube 9 is a negative control tube. Upon addition of the test organism suspension, final dilutions will be 1:20 through 1:2560 for tubes 1 through 8, respectively.
5. Add 0.5 mL of the test organism suspension to each of 9 tubes.
6. Add 0.5 mL of an appropriate **Difco** Listeria O Antigen to each of 9 tubes containing antiserum in the positive control row.
7. Shake the rack to mix. Incubate in a 50°C water bath for 2 h. Refrigerate overnight. Read for agglutination the following morning.

**User Quality Control:** At the time of use, test both positive and negative antigen controls to check performance of the antiserum, techniques and methodology. If results of the positive control or negative control are not as described, the test is invalid and results cannot be reported.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent NCCLS guidance and CLIA regulations for appropriate Quality Control practices.

#### RESULTS

##### Slide Test

1. Read and record results as follows:
 

4+	100% agglutination; background is clear to slightly hazy.
3+	75% agglutination; background is slightly cloudy.
2+	50% agglutination; background is moderately cloudy.
1+	25% agglutination; background is cloudy.
-	No agglutination
2. For the test isolate, a 3+ or greater agglutination within 1-2 min is a positive result.
3. A partial (less than 3+) or a delayed agglutination reaction should be considered negative.

##### Tube Test

1. Read and record results as follows:
 

4+	100% agglutination; background is clear to slightly hazy.
3+	75% agglutination; background is slightly cloudy.
2+	50% agglutination; background is moderately cloudy.
1+	25% agglutination; background is cloudy.
-	No agglutination
2. For the test isolate, a positive result is indicated by a 2+ or greater agglutination at 1:320.

#### LIMITATIONS OF THE PROCEDURE

1. Serological techniques employing **Difco** Listeria O Antisera serve as corroborative evidence for the identification of *L. monocytogenes*. Final identification cannot be made without consideration of morphological, serological and biochemical characterization.
2. Excessive heat from external sources (hot bacteriological loop, burner flame, light source, etc.) may prevent making a smooth suspension of the microorganism or cause evaporation or precipitation of the test mixture. False positive reactions may occur.
3. Rough culture isolates occur and will agglutinate spontaneously, causing agglutination of the negative control (autoagglutination). If autoagglutination occurs, the culture is rough and cannot be tested. Subculture to a noninhibitory medium, incubate and test the organism again.
4. Agglutination reactions of 3+ or greater in the slide test are interpreted as positive reactions. Cross-reactions resulting in a 1+ or 2+ agglutination are likely since there are somatic antigens shared among different organisms such as staphylococci, enterococci and *Bacillus* species.<sup>10</sup>
5. Prolonged exposure of reagents to temperatures other than those specified is detrimental to the products.
6. Exposure of Listeria O Antigens to temperatures below 2°C can result in autoagglutination. Antigens must be smooth uniform suspensions. Examine antigen vials for agglutination before use. Suspensions with agglutination are not usable and should be discarded.
7. It is important in this test to use the recommended time and temperature of incubation. Also, care should be taken to make certain that the water bath is in a location free of mechanical vibration.
8. Discard any **Difco** Listeria O Antiserum that is cloudy or has a precipitate after rehydration or storage.

**PERFORMANCE CHARACTERISTICS**

1. The positive control, using the appropriate *Listeria* O Antigen, should show 3+ or greater agglutination in the slide test.
2. The positive control, using the appropriate *Listeria* O Antigen, should show 2+ or greater agglutination at 1:320 in the tube test.
3. The negative control should show no agglutination.

Sensitivity of **Difco** *Listeria* O Antisera Types 1, 4 and Poly is determined by demonstrating appropriate reactivity in both the slide and tube tests, as defined in the "Results" section using the respective homologous **Difco** *Listeria* O Antigen. Specificity is determined by demonstrating non-reactivity against non-related (heterologous) **Difco** *Listeria* O Antigens.

**AVAILABILITY**

Cat. No.	Description
223001	<b>Difco</b> <sup>TM</sup> <i>Listeria</i> O Antiserum Type 1, 1 mL
223011	<b>Difco</b> <sup>TM</sup> <i>Listeria</i> O Antiserum Type 4, 1 mL
223021	<b>Difco</b> <sup>TM</sup> <i>Listeria</i> O Antiserum Poly, 1 mL
223031	<b>Difco</b> <sup>TM</sup> <i>Listeria</i> O Antigen Type 1 (Slide), 5 mL
223051	<b>Difco</b> <sup>TM</sup> <i>Listeria</i> O Antigen Type 1 (Tube), 25 mL
223041	<b>Difco</b> <sup>TM</sup> <i>Listeria</i> O Antigen Type 4 (Slide), 5 mL
223061	<b>Difco</b> <sup>TM</sup> <i>Listeria</i> O Antigen Type 4 (Tube), 25 mL
223143	<b>Difco</b> <sup>TM</sup> FA Buffer, Dried, 6 x 10 g
223142	<b>Difco</b> <sup>TM</sup> FA Buffer, Dried, 100 g
223141	<b>Difco</b> <sup>TM</sup> FA Buffer, Dried, 10 kg

**REFERENCES**

1. Murray, E.G.D., R.A. Webb and M.B.R. Swann. 1926. A disease of rabbits characterized by large mononuclear leucocytosis caused by a hitherto undescribed bacillus *Bacterium monocytogenes* (n. sp.). *J. Path. Bact.* 29:407-439.
2. Monk, J.D., R.S. Clavero, L.R. Beuchat, M.P. Doyle and R.E. Brackett. 1994. Irradiation inactivation of *Listeria monocytogenes* and *Staphylococcus aureus* in low and high fat, frozen and refrigerated ground beef. *J. Food Prot.* 57:969-974.
3. Wehr, H.M. 1987. *Listeria monocytogenes* – a current dilemma special report. *J. Assoc. Off. Anal. Chem.* 70:769-772.
4. Bremer, P.J., and C.M. Osborne. 1995. Thermal death times of *Listeria monocytogenes* in green shell mussels (*Perna canaliculus*) prepared for hot smoking. *J. Food Prot.* 58:604-608.
5. Patel, J.R., C.A. Hwang, L.R. Beuchat, M.P. Doyle and R.E. Brackett. 1995. Comparison of oxygen scavengers for their ability to enhance resuscitation of heat injured *Listeria monocytogenes*. *J. Food Prot.* 58:244-250.
6. Ryser, E.T., and C.W. Donnelly. 1992. *Listeria*, p. 343-356. In F.P. Downes and K. Ito (ed.), *Compendium of methods for the microbiological examination of foods*, 4th ed. American Public Health Association, Washington, D.C.
7. Seeliger, H.P.R., and K. Hohne. 1979. Serotyping of *Listeria monocytogenes* and related species, p. 31-49. In T. Bergen and J.R. Norris (ed.), *Methods in microbiology*, vol. 13. Academic Press, London, England.
8. Bille, J., J. Rocourt and B. Swaminathan. 2003. *Listeria* and *Erysipelothrix*, p. 461-471. In P.R. Murray, E.J. Baron, J.H. Tenover, M.A. Tenover and R.H. Tenover (ed.), *Manual of clinical microbiology*, 8th ed. American Society for Microbiology, Washington, D.C.
9. Pezzlo, M. 1994. Aerobic bacteriology, p. 1.0.1.-1.20.47. In H.D. Isenberg (ed.), *Clinical microbiology procedures handbook*, vol. 1. American Society for Microbiology, Washington, D.C.
10. Forbes, B.A., D.F. Sahn, A.S. Weissfeld and E.A. Trevino. 2002. *Bailey & Scott's diagnostic microbiology*, 11th ed. Mosby-Year Book, Inc., St. Louis, Mo.
11. Association of Official Analytical Chemists. 2001. *FDA Bacteriological analytical manual online*. <<http://www.cfsan.fda.gov/~ebam/bam-mm.html>>.
12. Flowers, R. S., W. Andrews, C. W. Donnelly, and E. Koenig. 1993. Pathogens in milk and milk products. In Marshall, R. T., (ed.), *Standard methods for the examination of dairy products*, 16th ed. American Public Health Association, Washington, D.C.

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© 2005 BD.
Difco is a trademark of Difco Laboratories, Inc., a subsidiary of Becton, Dickinson and Company.
BD and BD Logo are trademarks of Becton, Dickinson and Company.
BENEX Limited